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
**IMPLANT-retained suprastructures
passivated to a stress free fit through spark erosion**

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Passive fit with SAE Secotec Spark Erosion Technology

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See: **Quintessenz Journal of Dental Technology**

Volume 2 / Number 2 / April 2004, page 158-181

Edentulous and partially dentate patients: Quality assurance through interdisciplinary cooperation
Rübeling / Freesmeyer


JOMI, Vol. 19, No. 6/2004, page 810-818

Improving the fit of Implant-Supported Superstructures using the Spark Erosion Technique

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
* These system parts are used only for making dental prostheses!

	Seite/Page		Seite/Page
SECOTEC Informationen			
SECOTEC Information			
Sheffield-Test	4/5	Nobel Biocare UCLA - rotationsgeschützt	32/33
SAE-Funkenerosionsverfahren		Nobel Biocare UCLA-RP (Regular Platform)	34/35
<i>SAE Spark Erosion Method</i>	6/7	Nobel Biocare UCLA-NP (Narrow Platform)	34/35
Secotec-Verfahren unter Nutzung der Funkenerosion		Nobel Biocare UCLA-WP (Wide Platform)	34/35
<i>Secotec System using Spark Erosion</i>	8/9	Nobel Biocare UCLA-WP (Wide Platform)	34/35
SAE-Secotec Modellvorbereitung		PALTOP Advanced Multi Unit	34/35
<i>SAE Secotec Model Preparation</i>	10/11	Schütz Dental Group Schütz-Impla	36/37
Das formkonstante Epoxydharzmodell		Steri-Oss PME	36/37
<i>The consistently shaped epoxy resin model</i>	11-13	Straumann Bone Level Ø 3,5	36/37
Grafische Darstellung des Modellaufbaus		Straumann Bone Level Ø 4,5	36/37
<i>Diagram of Model Structure</i>	12/13	Straumann Bone Level Ø 6,5	38/39
SAE-Horizontalschleuder		Straumann Octa	38/39
SAE Horizontal Centrifuge	11/13	Straumann SynOcta	38/39
SAE-Präzisions-Parallelometer		Straumann SynOcta 1,5	38/39
SAE Precision Parallelometer	12/13	Straumann SynOcta WN (Wide Neck)	40/41
SAE-Funkenerosionsmaschine		Straumann Bone Level Ø 3,5	40/41
SAE Spark Erosion Machine	14/15	Verschraubtes Sekundärteil	
SAE-Grundausrüstung		Straumann Bone Level Ø 4,6	40/41
SAE Basic Equipment	16/17	Verschraubtes Sekundärteil	
		SPI Thommen Ø 3,5	40/41
SAE-SECOTEC-PROGRAMM 2017		SPI Thommen Ø 4,0	42/43
SAE-SECOTEC-PROGRAM 2017		SPI Thommen Ø 4,5	42/43
Ankylos	18/19	SPI Thommen Ø 5,0	42/43
Ankylos Balance Basisaufbau schmal	18/19	SPI Thommen Ø 6,0	42/43
Astra Tech Dental 20°	18/19	Zimmer Spline: Ø 3,25	44/45
Astra Tech Dental 45°	18/19	Zimmer Spline: Ø 4,0	44/45
Astra Tech Dental Uni Abutment	20/21	Zubehör SAE-Horizontalschleuder	
Bego-Semados	20/21	SAE-Präzisions-Parallelometer	44/45
Bego Multi Plus	20/21		
Biomet 3i Certain Ø 4,0 mm	20/21	SECOTEC-Hilfswerkzeug	
Biomet 3i Certain Ø 5,0 mm	22/23	SECOTEC Accessories	
Biomet 3i Certain Ø 6,0 mm	22/23	Drehmomentschlüssel und Rätsschrad	
Biomet 3i Standard Abutment AB 200-700	22/23	<i>Torque wrench and wrench wheel</i>	46/47
Biomet 3i Microminiplant (2.5 Hex)	22/23	Adapter für ISO-Schaft-Instrumente	
Biomet 3i Ø 6,0 UCLA WP (Wide Platform)	24/25	<i>Adapter for ISO shaft instruments</i>	46/47
Camlog/Conelog Ø 3,3 + Ø 3,8 + Ø 4,3	24/25	Konter- und Imbusschlüssel	
Camlog/Conelog Ø 5,0 + Ø 6,0	24/25	<i>Counter wrench and hex key</i>	46/47
Dentegris Ø 3,75	24/25	Gewindehülse und Schraube aus Titan	
Dentegris Ø 4,5+ Ø 4,75	26/27	<i>Titanium abutment and screw</i>	46/47
Friatec Frialit MP2Ø 3,4 + Ø 3,8 / Ø 4,5 / Ø 5,5	26/27		
MIS Multi Unit	26/27	Zur Patienteninformation	
MIS Multi Unit 2016	26/27	For Patient's Information	
MedentiBASE	28/29	Verschraubungen*	
Neoss-Multi Abutment	28/29	<i>Screw connections*</i>	48/49
NobelActive NP Ø 3,5	28/29	Titanröhrchen für Röntgen- und Bohrschablone	
NobelActive NP Ø 4,3	30/31	<i>Titanium cylinder for radiographic and drilling jig</i>	48/49
Nobel Biocare Mirus-Cone Ø 4,8 RP		Zubehör für das formkonstante	
(Regular Platform) = multi-unit	30/31	Funkenerosionsmodell	
Nobel Biocare Replace Select NP	30/31	<i>Accessories for the consistently designed</i>	
Nobel Biocare Replace Select RP	30/31	<i>spark erosion model</i>	48/49
Nobel Biocare Replace Select WP	32/33	Information für Implantatarbeiten	
Nobel Biocare Replace Select Ø 6,0 mm	32/33	<i>Information on implant-retained restorations</i>	50
Nobel Biocare Standard	32/33	Ncm-Einstellungen am Drehmomentschlüssel	
		<i>Ncm Adjustments to the torque wrench</i>	51

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Sheffield Test on CAD/CAM-milled bar structure

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

*One outer screw is fastened
before spark erosion passivation of the structure
no stress-free fit of the structure*

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*One center screw is fastened
before spark erosion passivation of the structure
no stress-free fit of the structure*

The SAE Spark Erosion Method for Achieving a Passive Fit of Implant Meso and Suprastructures

The conventional methods of casting framework in the dental laboratory will never lead to a passive fit of the sub and suprastructure situated on osseointegrated implants.


However, to comply with the demands of dentistry to achieve a passive fit of the implant related framework, many variables are considered. The cast suprastructure is segmented and reconnected, either by soldering or welding. This procedure leads to new inaccuracies, if only to an insignificant extent. But when the ceramic veneering is applied to the multispanded restoration, the framework shows a stress producing potential due to the shrinkage of the ceramic during firing. Again, we find an inaccurate fitting fixture.

These inaccuracies resulting from conventional prosthetic techniques can be corrected after completion of the suprastructure by applying the SAE Secotec Spark Erosion Method.

With the SAE Secotec Spark Erosion Technique, short circuit impulses in an electrical field burst out particles from the metallic surface of the suprastructure, starting from the areas that first come into contact with the lab implants and continuing until all contours of the framework are equally spark eroded, thus achieving uniform contact.

A special model is needed to provide electrical contact during the spark erosion procedure. The model shells furnish the opportunity to exchange the implant replicas for identical electrodes assuring an even flow of electricity and a sufficient production of sparks.

The Secotec System is applicable to all conductible metals and alloys and to acrylic or ceramic veneered suprastructures.

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**All screws are fastened
before spark erosion passivation of the structure
no stress-free fit of the structure**

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**All screws fastened
stress-free fit after spark erosion passivation
of the structure**

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What is Spark Erosion?

Spark Erosion is a method for shaping and eroding electrically conductive metals by controlled electrical current impulses between tool electrode (implant electrode) and work piece (implant meso or suprastructure) in the presence of a non-conductive (dielectric) fluid.

Compared with the machining procedure (milling – drilling), the erosion is an additional procedure for processing metals, however with the main difference that the rotation of the tool (drill cutter) is omitted and through the vertical propulsion of the fixed electrode, specifically for restorations on implants, it is possible to obtain exact fits with a customized design.

For 20 years spark erosion has been used in dental technology in the area of combined prosthodontics. With implant prosthodontics the procedure has been used for several years for passivation of the fit of meso and suprastructures on implants.


The contraction caused by the casting of implant meso and suprastructures can be quickly rectified by the spark erosion of the contact areas causing tension so that a passive fit of the meso or suprastructures on the implants can be obtained.

The principle of Spark Erosion

The work piece (meso or suprastructure) and electrode are put into the spark erosion machine and brought into the working position (illustration 1 and 2) so that they do not touch each other. A gap remains which is filled with an insulating liquid (dielectric fluid). Processing occurs in a container. The work piece and electrode are connected by an electrical cable to a direct current source. The power supply is interrupted by a switch. If this is closed, short circuit impulses occur between the work piece and electrode which give off flashovers from the electrode to the work piece.


These flashovers cause metal removal from the work piece (structure) on the areas hit first. Metal removal continues until all areas of the structure have full contact with the electrodes. The work piece sits fully circularly on the implant electrode whose surface design corresponds exactly to that of the implants or the abutments.

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
An electric field is developed between electrode and work piece by turning on the voltage.

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
The electrically conducted particles concentrate in the field between electrode and work piece and form a bridge.

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
Suddenly a channel of discharge develops over this bridge.

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
As discharge begins, the number of the electrically conducted particles increases like an avalanche, the electricity begins to flow.

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
Counter-pressure arises because of the viscosity of the dielectric fluid which ties up the channel of discharge between electrode and work piece. The concentrated heat brings the surface areas to melting.

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
The temperature increases, the melting continues into the surface areas. High temperatures are reached temporarily due to the vaporization of the melt.

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The flow of electricity is interrupted by switching off the voltage, the channel of discharge fails. This leads to sudden centrifugation and vaporization of the melt.

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There is fine electrode material in the dielectric fluid at the end of the discharge. The surface areas of the work piece and the electrode have each an ablation crater.

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Purpose of Dielectric Fluid

Insulation:


The dielectric fluid has the important function of insulating the work piece from the electrode. The electrical disruptive discharge must occur with as little a gap as possible. This increases the ablation capacity and the precision of the shape.

Cooling:

As the electric spark hits the work piece at about 10,000°C the dielectric fluid must cool the electrode and work piece. Overheating the electrode must be avoided so that it does not lead to augmented wear and tear. The metal vapours that occur during the spark erosion process must be able to condense in the liquid, and in addition, the resulting abla- tion particles must be washed away from the erosion zone to avoid any process errors.

SAE-SECOTEC – Das formkonstante SAE Modell

Modellvorbereitung – Step by Step / Model Preparation – Step by Step

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
Figs. 1 + 1a

The customized tray holding the impression in Impregum with the impression copings.

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Figs. 2 + 2a

The lab implant replicas of a matching implant system are attached to the Secotec model shells and then screwed to the impression copings embedded in the impression. All tightening procedures are performed with a torque wrench (item no. 82-0521) and a counter wrench (item no. 82-0531), hereby complying with the specific torque – Ncm. – Copper wire (item no. 82-0500) is attached to each model shell, thus connecting them to form an electrical circuit.

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Fig. 3

An effusion of elastic die acrylic is carried out. The implant replicas are completely covered with the soft acrylic and only the model shells are visible.

SAE-SECOTEC – SAE Model with consistent design

Modellvorbereitung – Step by Step / Model Preparation – Step by Step

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Fig. 5

The entire area is lightly dusted with silver powder to insulate against the Epoxy Resin.

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Fig. 6

A partial effusion is made with low contracting – 0.00 - 0.03 mm – SAE-Epoxy Resin (item no. 40-1060) and with the SAE Horizontal Centrifuge (item no. 55-0780) compression occurs.

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

Fig. 7

After all of the excess wax is removed the remaining effusion is completed in a type IV SAE Special Die Stone (item no. 70-1121 (light brown) / 70-1124 (ivory).

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Model Structure

- 10 die model
- 10a removable elastic portion
- 12 area for receiving implant replica and electrode
- 12a model shell
- 14 thread of model shell
- 16 copper wire contact
- 18 copper wire for conductivity (anode)



Abbildung 8

- a) Drehmomentschlüssel (82-0521) Einstellung 20 Ncm
- b) Eindrehwerkzeug
- c) Modellimplantat
- d) Modellhülse
- e) Konterschlüssel (82-0531)

Fig. 8

- a) torque wrench (82-0521) Setting 20 Ncm
- b) driver
- c) implant replica
- d) model shell
- e) counter wrench (82-0531)

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Basic Equipment for SAE Spark Erosion Machine

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Erosion Accessories

1. Chuck holder to receive chucks
Chucks to receive copper electrodes, chuck capacity 0.5-6.0 mm
2. Chuck holder 1.0-0,5 mm
3. Chuck holder 2.0-1.0 mm
4. Chuck holder 3.0-2.0 mm
5. Chuck holder 4.0-3.0 mm
6. Chuck holder 5.0-4.0 mm
7. Chuck holder 6.0-5.0 mm
8. Holder to engage superstructure for erosive machining with implant abutments
9. Universal electrode holder electrode holder for all T attachment erosions
10. Holder for centering system to achieve exact 0 position in machine
11. Holder for all latch attachment erosions
12. Depth gauge to enable the adjustment of original depth for finishing erosive process
13. Electrode guide for latch axis
14. Angle holder latch axis Dielectric fluid (not shown) Sleeves for electrode guide for 30-1602 (not shown)

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Accessories

15. Model mount for EDM 2005 and EDM 2000
16. Angled plate for model mount
17. Securing clamps for model (1 pair)
18. Clamp for model Rubber circulars 0.5 - 2.5 mm for EDM 2005 and EDM 2000 (not shown)

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
► Ankylos

1. Model shell standard
2. Implant replica
3. Implant electrode 4,45 Ø
4. Plastic burn-out cylinder
5. Prosthetic screw ❶
6. Waxing screw
7. Impression screw
8. Imprint cap
9. Driver implant replica
10. Driver electrode and screw
11. Driver waxing screw

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► Ankylos Balance

1. Model shell standard
2. Implant replica
3. Implant electrode
4. Plastic burn-out cylinder
5. Prosthetic screw ❶
6. Impression screw
7. Imprint cap
8. Driver implant replica
9. Driver electrode and screw

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► Astra 45°

1. Model shell standard
2. Implant replica
3. Implant electrode
4. Plastic burn-out cylinder
5. Prosthetic screw ❶
6. Driver implant replica
7. Driver electrode and screw

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Astra Abutment

1. Model shell standard
2. Implant replica
3. Implant electrode
4. Plastic burn-out cylinder
5. Prosthetic screw ①
6. Driver electrode, implant and screw

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Bego-Semados

1. Model shell standard
2. Implant replica
3. Implant electrode
4. Plastic burn-out cylinder
5. Prosthetic screw ①
6. Waxing screw
7. Driver implant replica
8. Driver electrode and screw

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3i Certain Ø 4.0 mm

1. Model shell standard
2. Implant replica Ø 4.0 mm
3. Implant electrode Ø 4.0 mm
4. Plastic burn-out cylinder Ø 4.0 mm
5. Prosthetic screw hexagon ①
6. Prosthetic screw square ①
7. Driver implant replica
8. Driver electrode and screw
9. Driver screw square

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3i Certain Ø 5.0 mm

1. Model shell standard
2. Implant replica Ø 5.0 mm
3. Implant electrode Ø 5.0 mm
4. Plastic burn-out cylinder Ø 5.0 mm
5. Prosthetic screw hexagon ❶
6. Prosthetic screw square ❶
7. Driver implant replica
8. Driver electrode and screw
9. Driver screw square

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3i Certain Ø 6.0 mm

1. Model shell standard
2. Implant replica Ø 6.0 mm
3. Implant electrode Ø 6.0 mm
4. Plastic burn-out cylinder Ø 6.0 mm
5. Prosthetic screw hexagon ❶
6. Prosthetic screw square ❶
7. Driver implant replica
8. Driver electrode and screw
9. Driver screw square

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► Nobel Biocare Replace Select WP

1. Model shell standard
2. Implant replica
3. Implant electrode
4. Plastic burn-out cylinder
5. Prosthetic screw 2,0 ①
6. Driver implant replica
7. Driver electrode and screw

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► Nobel Biocare Standard

1. Model shell standard
2. Implant replica
3. Implant electrode
4. Plastic burn-out cylinder
5. Prosthetic screw ①
6. Driver implant replica
7. Driver electrode and screw

This system is being phased out and will be replaced by Nobel Biocare Multi Unit. Identical with Mirus-Cone Ø 4,8 RP.

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► Nobel Biocare UCLA

1. Model shell
2. Implant replica
3. Implant electrode
4. Plastic burn-out cylinder
5. Prosthetic screw ①
6. Direction screw
7. Driver electrode and screw

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► Straumann SynOcta

1. Model shell standard
2. Implant replica
- * SynOcta secondary part
(original system component/Abutment)
3. Implant electrode
4. Plastic burn-out cylinder
5. Prosthetic screw ❶
6. Driver implant replica
7. Driver SCS screw
8. Driver electrode and screw

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► Straumann SynOcta 1.5

1. Model shell standard
2. Implant replica
- * SynOcta 1.5 secondary part
(original system component/Abutment)
3. Implant electrode
4. Plastic burn-out cylinder
5. Prosthetic screw ❶
6. Driver implant replica
7. Driver SCS screw
8. Driver electrode and screw

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SAE-SECOTEC-Implantatsystem

SAE-SECOTEC Sistema Implantare

SAE-Secotec für **Straumann SynOcta WN** (Wide Neck)

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
► Straumann SynOcta WN

1. Model shell standard
2. Implant replica
- * SynOcta WN secondary part
(original system component/Abutment)
3. Implant electrode
4. Plastic burn-out cylinder
5. Prosthetic screw ①
6. Driver implant replica
7. Driver SCS screw
8. Driver electrode and screw

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► Straumann Bone Level Ø 3,5

1. Model shell standard
2. Implant replica*
3. Implant electrode 3.5 Ø
4. Plastic burn-out cylinder
5. Prosthetic screw ①
6. Driver implant replica
7. Driver electrode and screw

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► Straumann Bone Level Ø 4,6

1. Model shell standard
2. Implant replica
3. Implant electrode 4,6 Ø
4. Plastic burn-out cylinder
5. Prosthetic screw ①
6. Driver implant replica
7. Driver electrode and screw

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► SPI Thommen Ø 3.5

1. Model shell standard
2. Implant replica
3. Implant electrode
4. Plastic burn-out cylinder
5. Prosthetic screw M 1.4 ①
6. Driver implant replica
7. Driver electrode and screw

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► SPI Thommen ø 4,0

1. Model shell standard
2. Implant replica
3. Implant electrode
4. Plastic burn-out cylinder
5. Prosthetic screw M 1.6 ①
6. Driver implant replica
7. Driver electrode and screw

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Implant-related screws of the Secotec System manufactured by SAE are for laboratory use only. To clinically secure implant-retained mesostructures and superstructures please use the original system screws – no SAE Secotec screws – risk of breakage!

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► Zimmer Spline: ø 3.25

1. Model shell
2. Implant electrode GH 1
3. Implant electrode GH 2
4. Implant electrode GH 3
5. Plastic burn-out cylinder
6. Prosthetic screw ❶
7. Driver electrode and screw

* original system component/Abutment

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Torque wrench and Wrench wheel

1. SAE-Wrench wheel for Camlog®-SAE-Torque wrench
2. Camlog®-Wrench wheel for Camlog®-SAE-Torque wrench
3. Camlog®-SAE-Torque wrench

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Adapter for ISO shaft instruments

1. Adapter for ISO shaft instruments

The adapter gives you the opportunity to insert any tool with ISO-shaft into the torque wrench.

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Counter wrench and hex key

1. Counter wrench for model shell standard
2. Hex key for model shell
3. Camlog®-SAE-Torque wrench

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Titanium abutment and screw

1. Cu electrode for first erosion Ø 1,5 for screw
2. Electrode for titanium first erosion Ø 1,5 for screw
3. Cu electrode for second erosion Ø 2,9 for abutment
4. Titanium abutment
5. Transverse tight screw titanium
6. Abutment Au
7. Transverse tight screw Au
8. Driver

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Screw Connections

1. Cylinder head screw titanium for threaded sleeve, small, M 1.6 mm
2. Grub screw titanium for threaded sleeve, small, M 1.6 mm
3. Threaded sleeve titanium, small, M 1.6 mm
4. Threaded sleeve CoCrMo, small, M 1.6 mm
5. Cylinder head screw large, titanium, M 2.0 mm
6. Grub screw titanium for titanium threaded sleeve, large, M 2.0 mm
7. Grub screw gold for gold threaded sleeve, large M 2.0 mm
8. Threaded sleeve CrCoMo, large, M 2.0 mm
9. Threaded sleeve titanium, large, M 2.0 mm
10. Threaded sleeve gold, large, M 2.0 mm
11. Driver
12. Driver - hand driver

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*These system parts are used only for making dental prostheses!

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Titanium cylinder

1. Titanium cylinder inner Ø 3.1 for pilot drill
2. Titanium cylinder inner Ø 2.00 for Lindemann miller length 8 mm, fits into 82-0497
3. Titanium cylinder inner Ø 2,4

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Mixing Ratio 1 to 6:

volume: hardener 1 ml: Resin 6 ml
weight: hardener 1 g: Resin 10 g
(see also QZ 3/2004, pages 238 - 250)

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Information on implant-retained restorations

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Information on implant-retained restorations – concerning screws ①

Implant-related screws of the Secotec System manufactured by SAE are for laboratory use only. To clinically secure implant-retained mesostructures and superstructures please use the original systems screws – no SAE Secotec screws – risk of breakage!

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All SAE-Secotec accessories are for laboratory use only.

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
Ncm Adjustments to the Torque Wrench

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

Avoiding mishaps with implant-retained prostheses.
Passive fit using spark erosion!

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